

TYPES OF RAMPS

Twin monoramp

Lenght [cm]	Weight [kg]	Max. load [kg]
RP 100	4,1	650
RP 120	5,1	540
RP 150	6,2	430
RP 200	8,3	320



Ideal for uniform height differences, where there is a need for the least possible weight, and where the ramp is used often in the same place

Twin telescopic ramp

Lenght [cm]	Weight [kg]	Max. load [kg]
RZ 100/65	5,0	650
RZ 125/78	6,2	540
RZ 150/90	8,5	430
RZ 200/120	10	320



A flexible ramp that can be used in many places and with varying height differences where height and lenght are variable. Takes up little space when not in use.

Monoramp


Lenght [cm]	Weight [kg]	Max. load [kg]
RSP 100	9,7	325
RSP 125	12,1	325
RSP 150	14,4	325
RSP 200	19,1	325
RSP 300	28,6	200



Ideal in places where the ramp is permanently in position or always used in the same place. Provides increased safety for users and helpers. Good for scooters and electric wheelchairs where comfort is a major consideration.

Folding monoramp

Lenght [cm]	Weight [kg]	Max. load [kg]
RSZ 150/78	17,5	325
RSZ 200/103	22,2	325
RSZ 250/128	26,9	200
RSZ 300/153	31,7	200



Ideal when the ramp is used in the same place and for placing in the car. Provides increased safety for users and helpers. Good for scooters and electric wheelchairs where comfort is a major consideration. Takes up less space than plain ramps when not in use.

EasyFold ramp


Lenght [cm]	Weight [kg]	Max. load [kg]
RSS 90	7,3	325
RSS 120	9,6	325
RSS 150	12,1	325
RSS 180	14,3	325



Suitable for minor height differences. Folds along its width to half its original size and provides a wide surface in its unfolded state. Can be carried and transported as easily as a suitcase. Excellent for scooters and electric wheelchairs.

Telescopic EasyFold ramp

Lenght [cm]	Weight [kg]	Max. load [kg]
RZS 180/112	16,4	325
RZS 210/127	18,9	325
RZS 240/143	21,4	325
RZS 300/173	26,4	200



Perfect for overcoming larger height differences. A wide ramp that is easily extended to the desired lenght. Can be folded along its width and extended in lenght, thus taking up less space than other long ramps - also during transport and storage. Excellent for scooters and electric wheelchairs.